

MS® SimplePak™

Gridded Membrane Filter

Sterile, continuous Packaging,

dispensed with contamination-free

for microbiological analysis

MS® SimplePak™ are available pore size ranging from 0.2 μ m to 0.8 μ m, which are made from mixed esters of cellulose; with two colors: white and black.; two diameters in either: 47mm and 50mm(0.45 μ m and 0.80 μ m).

- 1 Excellent retention and colony growth, high recovery rates of microorganisms
- 2 Sterile by gamma irradiation
- 3 High-throughput
- 4 Traceability: filter type, diameter, lot number, and expiration date are all printed on each box
- 5 “Multi-fit” are designed as “one size fits all”, which can be used in all standard dispensers.



Wide Range of Application

MS® SimplePak™ can be used in general purpose of all microorganisms examination in medical, wine, beverage and pharmaceutical industry, designed for use in vacuum and pressure filtration applications.

White membrane with black grid

Applications:

Analysis of E. Coli bacteria

Examination of:


Dark colonies


Pore Size:

0.22 μ m, 0.45 μ m, 0.8 μ m

Diameter:

47mm, 50mm





Black membrane with white grid:

Applications:

Analysis of yeast and mold and Legionella sp.

Examination of:

White, being colonies

Pore Size:

0.45µm, 0.8µm

Diameter:

47mm, 50mm

The easiest way to work with MS® Simple Pak™

Convenient dispenser

- 1 More easier to open filter packages
- 2 Save time
- 3 Reduce the risk of membrane contamination
- 4 Easily moved to everywhere you want



In the process of using dispenser, you can only press the dispenser's lever, the sterile membrane is unwrapped from its package, and your figure never need to close to the filter. The dispenser will make your work more efficient and more safe!

Ordering information

MS® Simple Pak™ MCE gridded membrane filter

Cat No.	Filter Color	Pore Size (µm)	Filter Diameter (mm)	Pcs/pk
MFMCCE047022CW	White with black grid	0.22	47	150
MFMCCE047045CW		0.45	47	150
MFMCCE047080CW		0.80	47	150
MFMCCE050022CW		0.22	50	150
MFMCCE050045CW		0.45	50	150
MFMCCE050080CW		0.80	50	150
MFMCCE047045CB	Black with white grid	0.45	47	150
MFMCCE047080CB		0.80	47	150
MFMCCE050045CB		0.45	50	150
MFMCCE050080CB		0.80	50	150

Membrane Dispenser

Description	Pcs/pk	Cat No.
Membrane Dispenser	1 each	MFDISP001

Membrane Solutions LLC

USA

Toll free: 866-528-4572

Fax: 732-412-4040


Address: 1108 Spring View Lane Plano
TX75075

China

Toll free: 400-620-1660

Fax: 0086-21-51687551

Address: 2201, No.1759 North
Zhongshan Road, Shanghai



For more enquiry:

Please write to: info@membrane-solutions.com